



## Spirometer disposables

### High quality disposables at affordable prices.

MD Spiro offers a complete line of spirometer disposables designed to meet or exceed all manufacturers' specifications.

All MD Spiro disposables are manufactured from the highest quality materials, packaged in convenient quantities and offered at the best prices available today.

MD Spiro also offers a complete line of filter adapters that fit every competitor's spirometer.

From the low cost standard cardboard mouthpiece to the high efficiency, yet affordable SpiroSafe filter, MD Spiro has the spirometer disposables you need.

Call **1-800-588-3381** to place an order and cut your disposable costs today.



#### SPIROSAFE™ FILTER #3385

##### Features:

- Virus Filtration Efficiency 99.84%
- Bacterial Filtration Efficiency 99.99%
- Flow Resistance .48 cm/H<sub>2</sub>O/L/Sec.
- Dead Space 55 cc
- Electrostat™ Filter Media
- 30 mm OD x 27 mm ID
- 100/bag, 3 bags/case



#### MICROCHECK™ MOUTHPIECE #3395 (Cardboard) #3397 (Plastic)

##### Features:

- One Way Valve
- Peace of Mind for Patients
- Great for Peak Flow Meters
- Lower Air Resistance
- Expiratory Only Testing
- 30 mm OD
- 200/box



#### DISPOSABLES # various

##### Features:

- High Quality Products
- Lowest Prices Available
- Fit All Competitors' Units
- Convenient Package Quantities

#### MD Spiro

803 Webster Street  
Lewiston, ME 04240  
Telephone 800-588-3381  
Fax 207-786-7280  
[www.mdspiro.com](http://www.mdspiro.com)