

## Viral/Bacterial Mouthpiece for CO Breath Testing

Micro Direct is pleased to introduce the **NEW SafeBreath™** viral/bacterial filtered mouthpiece for carbon monoxide (CO) breath testing.

**SafeBreath™** is a unique solution offering a combined filter and cardboard mouthpiece with >99% protection against both bacteria and viruses.

Using **SafeBreath™** not only protects your patient from potential cross infection during testing but also protects and increases the life of your fuel cell and eliminates the need for and cleaning of the plastic mouthpiece adapter!

The **SafeBreath™** (order #FM200) is packaged in a handy dispenser box for easy storage and dispensing.



### Features and Specifications

- Disposable
- Cleaner and safer
- >99% protection against both bacteria and viruses
- Eliminates plastic mouthpiece adapters saving money
- No more cleaning of plastic adapters saving time
- 22 mm OD
- Order #FM200 (\$70/box of 200)