

Affordable Spirometer For Home Use

The Micro Direct **spiro**✓home spirometer was specifically designed for situations where low cost, easy to use, accurate spirometry is required.

The **spiro**✓accurately measures, displays and stores the actual FEV₁, FEV₆ and the percent of Personal Best FEV₁. The personal best FEV₁ can be set by the clinician with future patient efforts compared to their personal best. The three color action zones are preset but can be customized by the clinician.

Features

- Measures FEV₁, FEV₆ and % of Personal Best FEV₁
- Simple to use
- Large, easy to read display
- Quality blow indicator
- FEV₁ zones can be personalized
- Stores 200 test sessions
- Easy to clean, detachable flow head
- Light, compact, attractive; makes it convenient to carry



Specifications

Product:

Spiro✓Home Monitor

Parameters Displayed:

FEV₁, FEV₆ and % Personal Best FEV₁

Memory:

200 Test Sessions

Accuracy:

Better than +/- 3% to ATS/ERS Standards

Volume Range:

0 - 9.99 Liters/Second (L/S) BTPS

Sensor:

Stator/Rotor

Flow Impedance:

Better than 0.15kPa/L/S at 14 L/S

Display:

Liquid Crystal Display (LCD)

Power Supply:

2 AAA Alkaline Batteries

Size:

4.5" x 2.5" x 2"

Weight:

Unit Only: 2.5 oz

Packed: 4.5 oz



Micro Direct, Inc.

803 Webster Street
Lewiston, ME 04240

Telephone 800-588-3381

Fax 207-786-7280

Direct www.mdspiro.com