

NLHEP QC Grades

The quality of each session is graded according to the following criteria:

QC Grade	Criteria
A	At least two acceptable maneuvers, with the largest two FEV ₁ 's matching better than 100 ml and the largest two FEV ₆ 's matching better than 100 ml
B	At least two acceptable maneuvers, with the largest two FEV ₁ 's matching better than 150 ml
C	At least two acceptable maneuvers, with the largest two FEV ₁ 's matching between 200 ml and 150 ml
D	Only one acceptable maneuver or two with the largest two FEV ₁ 's matching less than 200 ml
F	No acceptable maneuvers