

It's Important - Your office sees patients with indications for spirometry on a daily basis. Spirometry is the gold standard to support or exclude the diagnosis of both asthma and COPD and then for management once a diagnosis is made. The earlier you detect the disease; the easier you can deal with it.

It's Good Medicine - Your patients are correctly diagnosed and treated more accurately and conveniently in your office.

It's Easy - Unless the equipment is easy to use, it will not be fully utilized by your staff. All Micro Direct spirometers are menu driven via a large graphic display and the test takes less than five minutes to perform. The easy-to-read reports, quality checks and built-in interpretation assist with your diagnosis.

It's Safe - The SpiroSafe Filter is a disposable viral/bacterial filtered mouthpiece that helps to protect your patients and staff. It is >99% effective against viral and bacterial cross-contamination.

It's Good Business - Office spirometry is third party reimbursable with an average reimbursement of \$28 for a simple five minute non-invasive test. The average primary care physician tests ten patients per week and pays for their equipment in three months or less then generates over \$14,000 per year in additional revenue.

Why Micro Direct? - Only Micro Direct offers a complete line of automated spirometers designed to meet your needs. Each utilize inexpensive, readily available mouthpieces and meet all ATS/ERS standards. Micro Direct's spirometry experts will help you choose the right spirometer and then fully support your staff from proper testing techniques through billing questions; and all Micro Direct spirometers are backed by a 30-day money back guarantee.



SpiroSafe Filter
Order #3385
\$110/box of 100




MicroCheck One-Way Mouthpiece
Order #3395
\$80/box of 200



Adult Cardboard Mouthpiece
Order #3314SB or #3314B5
\$35/box of 200
\$75/box of 500



Pediatric Mouthpiece
Order #3301
\$17/bag of 100



asmaPLAN+ Peak Flow Meter
Order #43832
\$11 each
\$200/case of 25




Ask about discounted pricing on blanket orders for mouthpieces or filters

Spirometry Solutions from Micro Direct

Your Spirometry Specialist

<p>Alpha (MD6000) Desktop</p> <ul style="list-style-type: none"> ❖ Accurate, Stable Fleisch Pneumotach ❖ Updated Test Quality Grading System ❖ Meets 2019 ATS/ERS Standards ❖ Live Testing Prompts ❖ GLI Predicted Equations with LLN & %Pred ❖ Automatic Interpretation based on ATS/ERS & GOLD Standards ❖ Clear, Easy-to-Read Report ❖ Pre/Post Bronchodilator Comparison ❖ PDF Report Generation Software ❖ Autoclavable Flow Head ❖ Complete with Durable Carry Case ❖ Comes with Connect Software to provide connectivity to your EMR via HL-7 and GDT ❖ 2-Year or 5-Year Warranty if registered within 30-days of purchase <p>\$2,499</p>	<p>Pneumotrac (MD6800) PC-Based</p> <ul style="list-style-type: none"> ❖ Accurate, Stable Fleisch Pneumotach ❖ Supplied with SpiroTrac 6 PC Software ❖ Meets 2019 ATS/ERS Standards ❖ Real Time FV Loop or VT Curve ❖ Choice of Child Friendly Incentives ❖ 8 1/2" x 11" Printout ❖ Pre/Post Bronchodilator Comparison ❖ Automatic Interpretation based on ATS/ERS & GOLD Standards ❖ Clear, Easy-To-Read Reports ❖ GLI predicted equations and Z-scores. ❖ On-Screen Test Quality Prompts ❖ Autoclavable Flow Head ❖ Comes with Connect Software to provide connectivity to your EMR via HL-7 and GDT ❖ 2-Year or 5-Year Warranty if registered within 30-days of purchase <p>\$2,149</p>	<p>SpiroBank II (MD10) Stand Alone or PC Based</p> <ul style="list-style-type: none"> ❖ Measures and Displays FVC, FEV₁, FEV₁%, PEF, FEF₂₅₋₇₅, FET, ERV, VC, IVC and IC ❖ Choice of Predicted Values ❖ Pre/Post Bronchodilator Comparison ❖ Meets ATS/ERS 2019 Standards ❖ 10,000 Test Memory ❖ Long Lasting Lithium Rechargeable Battery via USB ❖ High Resolution Backlit Display ❖ Choice of Reusable or Disposable Turbine ❖ EHR/EMR Connectivity ❖ Comes with winspiroPRO PC Software ❖ 1-Year Parts and Labor Warranty <p>\$840</p>	<p>Spiro✓ Home (MD01) Home Screening</p> <ul style="list-style-type: none"> ❖ Measures and Displays FEV₁, FEV₆ and % of Personal Best FEV₁ ❖ Quick and Easy to Use ❖ Quality Blow Indicator ❖ Large Graphical Display ❖ Stores 200 Test ❖ 1-Year Warranty <p>\$295</p>
<p>Spiro✓ (MD02) Screening</p> <ul style="list-style-type: none"> ❖ Measures and Displays FEV₁, and FEV₁ % Predicted ❖ 'Lung Age' Estimation ❖ Quick and Easy to Use ❖ Large Graphical Display ❖ Standard 30 mm OD Mouthpiece ❖ 1-Year Warranty <p>\$175</p>			



Micro Direct, Inc. • 803 Webster Street, Lewiston ME 04240

Telephone 207-786-7808 • 1-800-588-3381 • Fax 207-786-7280 • www.mdspiro.com • e-mail: sales@mdspiro.com